

# ACTIVA 700

## 4 CHANNELS ELECTROTHERAPY | 271 PROGRAMS

PRODUCT FOR BEAUTY AND FITNESS, DESIGNED WITH SPECIFIC PROGRAMS FOR FIRING, TONING AND REJUVENATION. THE SPECIAL HANDPIECE G-TRODE (SUPPLIED) MAKES ACTIVA 700 A TOP OF THE RANGE OF THE ELECTROSTIMULATION LINE.

60 BEAUTY | 58 FITNESS | 53 SPORT | 36 SERIAL SEQUENTIAL STIMULATION | 21 SPECIALS SPORT | 13 SKIN-FACE  
12 MICROCURRENTS | 12 TENS | 3 REHAB | 3 INCONTINENCE

### TECHNICAL FEATURES

#### Display

Visible area size 2,6"

#### Channels

4 independent (8 electrodes)

#### Frequency

0,3-150Hz \*

#### Pulse amplitude

40-450µs \*

#### Power

0-120mA per channel

#### Power supply

Rechargeable batteries

#### Size

mm 100x160x35


#### Weight

gr 454

\* According to the type of current

### EQUIPMENT

- 1 Bag
- 1 Stimulator ACTIVA 700
- 1 Handpiece G-Trode
- 4 Cables for electrode connection
- 2 Cables for microcurrents and ionophoresis
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 1 Gel flask
- 1 Charger
- 1 Operating manual



CE  
0476



# ACTIVA 700

## 4 CHANNELS ELECTROTHERAPY | 271 PROGRAMS

### PROGRAMS LIST

#### FITNESS

Firming

Bio-Pulse firming

Sculpting

Bio-Pulse sculpting

Toning

Mass building

Body sculpting

Definition

Jogging

Anaerobic fitness

Aerobic fitness

Cramp prevention

#### BEAUTY

Drainage

Bio-Pulse drainage

Lipolysis

Post-pregnancy lipolysis

Bio-Pulse relaxing massage

Toning massage

Energizing massage

Connective massage

Swollen arms

Face capillaries

Definition

Slim figure

Skin tone improvement

Post-pregnancy drainage

Post-pregnancy firming

Breast firming

Breast sculpting

#### SPORT

Capillarization

Warm-up

Pre-competition warm-up

Active recovery

Maximum strength

Endurance strength

Explosive strength

Reactivity

Post-competition recovery

Hypertrophy

#### SPECIAL SPORT

Running

Skiing

Golf

#### REHAB

Quadriceps atrophy (with knee prosthesis)

Recovery after ACL surgery

Shoulder subluxation prevention

#### TENS and PAIN

Conventional antalgic tens

Modulated antalgic tens

Endorphinic tens

Menstrual pain

Cronich pain

Cervical pain

Muscle pain

Knee pain

Scapulohumeral periarthritis

Chronic lumbago

Bursitis-tendinitis

Osteoarthritis

#### G-PULSE

Breast microlifting

Cleavage microlifting

Face microlifting

Skin elasticity

Cellulite

BioSkin collagen

Wrinkles

Eye wrinkles

Expression wrinkles

Stretch marks

Nail strenghtening

Bio peeling

Hematoma

#### MICROCURRENTS

Epicondylitis

Scapulohumeral periarthritis

Contusion

Sciatica

Lumbago

Acute pain

Articular pain

Stiff neck

Whiplash

Shoulder sprain

Knee sprain

Tendon inflammation

#### INCONTINENCE

Mixed incontinence

Stress incontinence

URGE incontinence

#### SERIAL SEQUENTIAL STIMULATION

The "3S" program list includes 36 parameter combinations.

The "3S" programs are characterized by a delay activation of the channels 3 and 4 compared with the channels 1 and 2.

The Serial Sequential Stimulation permits to stimulate the musculature in kinetic chain thanks to the differentiated activation times of the muscular groups involved.